


DID YOU KNOW?

YOU CAN RECYCLE YOUR
MilliporeSigma CARTRIDGES?

Recycle Your Used MilliporeSigma Water System Cartridges

1. Drain used cartridge in sink for 24 hours
2. Deposit drained cartridge in the drum located in the **Biolabs 1041 "Reuse Room"** or request pickup as you would for chemical waste
3. Cartridges will be sent to a recycling facility and made into plastic lumber!

For more information please contact:
anthony_michetti@harvard.edu



In partnership with



HARVARD
Office for Sustainability



TRIUMVIRATE
ENVIRONMENTAL